




Livarian Labs

Research Reference Sheet

CJC-1295 (No Dac) + Ipa 5mg/5mg

Livarian Labs Research Reference Sheet CJC-1295 (No Dac) + Ipa 5mg/5mg • Reconstitution Calculation Reference

Mathematical conversion examples only. This sheet is for laboratory research reference and does not provide dosing, administration, therapeutic, diagnostic or clinical guidance.

Standard Reconstitution Example

Product	JC-1295 (No Dac) + Ipamorelin 5mg/5mg
Form	Lyophilised research compound
Compound amount	10.000 mg
Calculated concentration	3.333 mg/mL
Research-use position	Laboratory research use only
Calculator page	Research Calculator Hub
COA / Verification	Batch Verification / COA page
Diluent volume	3.000 mL
Approx. per 0.1mL	333.333 mcg
Document type	Reference sheet

Calculation Examples — 10mg vial + 3mL diluent

Vial amount	Diluent volume	Concentration
10 mg	3 mL	3.333 mg/mL
10 mg	3 mL	3.333 mg/mL
10 mg	3 mL	3.333 mg/mL
10 mg	3 mL	3.333 mg/mL

Target research amount	Calculated volume	U-100 syringe reference
100 mcg	0.030 mL	3.0 units
150 mcg	0.045 mL	4.5 units
200 mcg	0.060 mL	6.0 units
250 mcg	0.075 mL	7.5 units

Formula Reference

Concentration	compound amount ÷ diluent volume
Calculated volume	target research amount ÷ concentration
U-100 syringe reference	calculated mL × 100

Research Use Only

This reference sheet provides mathematical conversion examples only. It is not dosing advice and does not provide administration, clinical, therapeutic or diagnostic guidance. Livarian Labs products are supplied strictly for laboratory research use only. Not for human consumption.

Livarian Support Links

View Batch Verification / COA	/pages/batch-verification-coa
Use Research Calculator	/pages/calculator-research-info
Shop Research Compounds	/collections/all