



# Livarian Labs

Research Reference Sheet

TESAMORELIN 10mg

## Livarian Labs Research Reference Sheet TESAMORELIN 10mg • Reconstitution Calculation Reference

*Mathematical conversion examples only. This sheet is for laboratory research reference and does not provide dosing, administration, therapeutic, diagnostic or clinical guidance.*

### Standard Reconstitution Example

Product	TESAMORELIN 10mg
Form	Lyophilised research compound
Compound amount	10.000 mg
Calculated concentration	3.333 mg/mL
Research-use position	Laboratory research use only
Calculator page	Research Calculator Hub
COA / Verification	Batch Verification / COA page
Diluent volume	3.000 mL
Approx. per 0.1mL	333.333 mcg
Document type	Reference sheet

### Calculation Examples — 10mg vial + 3mL diluent

Vial amount	Diluent volume	Concentration
10 mg	3 mL	3.333 mg/mL
10 mg	3 mL	3.333 mg/mL
10 mg	3 mL	3.333 mg/mL
10 mg	3 mL	3.333 mg/mL

Target research amount	Calculated volume	U-100 syringe reference
500 mcg (0.5 mg)	0.15 mL	15 Units
1000 mcg (1.0 mg)	0.30 mL	30 Units
1500 mcg (1.5mg)	0.45 mL	45 Units
2000 mcg (2.0 mg)	0.60 mL	60 Units

### Formula Reference

Concentration	compound amount ÷ diluent volume
Calculated volume	target research amount ÷ concentration
U-100 syringe reference	calculated mL × 100

### Research Use Only

This reference sheet provides mathematical conversion examples only. It is not dosing advice and does not provide administration, clinical, therapeutic or diagnostic guidance. Livarian Labs products are supplied strictly for laboratory research use only. Not for human consumption.

### Livarian Support Links

View Batch Verification / COA	<a href="/pages/batch-verification-coa">/pages/batch-verification-coa</a>
Use Research Calculator	<a href="/pages/calculator-research-info">/pages/calculator-research-info</a>
Shop Research Compounds	<a href="/collections/all">/collections/all</a>